

KAT-PB-030

benefit  
quattro



**Instruction for Use**  
**BENEFIT® and quattro® systems**

E

E

REF

Warning

Caution

PSM

Recommendation

Temporary

Recommendation

3. INTENDED USE

4. INDICATIONS FOR USE

Distalization

5. INTENDED USER AND USE ENVIRONMENT

Warning

Caution

US Federal law

6. PATIENT POPULATION

7. CONTRAINDICATIONS

8. SAFETY PRECAUTIONS

9. UNDESIRABLE SIDE-EFFECTS AND COMPLICATIONS

- insufficient bone quantity or quality at the insertion site
- overloading
- debilitation
- use of a screw of inadequate diameter or length
- root contact during insertion
- manipulation by the patient's fingers or tongue
- poor oral hygiene
- non-compliance with orthodontic treatment
- application of excessive force on the appliance
- too long a lever arm, if the miniscrew is inserted in a region where the gingiva or mucosa is too thick
- peri-implantitis from insertion in the mucosal region
- insufficient primary stability
- orthodontist's inexperience regarding the treatment modality
- bone damage on insertion from stress or overheating
- bone density and cortical bone thickness

Furthermore, implant breakage could be caused by an excessive resistance during insertion or an inappropriate device selection.



Fig. 1: Bending precaution

**10. SURGICAL PROCEDURES**

Prior to surgery, evaluate the preferred implant site carefully for bone quality and quantity, using the lateral and anteroposterior cephalometric films, panoramic or periapical X-rays or computed tomographic scans. Select the insertion area and the screw size. Use local anesthesia. Drill a pilot hole through the cortical bone if applicable. Use the appropriate screwdriver (see accessories) to pick up the screw from the sterile holder and insert it directly through the gingiva. Screw insertion is done manually with steady rotational movements with the selected accessories until the screw heads touch the gingiva.

If a torque limited screwdriver is used, apply the following adjustment:

- 1.5 mm Screw max. 15 Nm,
- 2.0 mm Screw max. 30 Nm,
- 2.3 mm Screw max. 40 Nm.

The screw can be loaded immediately. Connect the screw with the selected orthodontic appliance. Screw removal is done by turning the screw counter-clockwise without local anesthesia. Following the insertion of the screws, the BENEFITs are bent either directly chairside or by means of impression taking and a plaster model in a laboratory. After precise bending of the BENEFITs to fit the implantable screw the wire is bent to the desired position (see safety precautions) according to the anatomical situation of the patient without touching the soft tissue. If the BENEFIT is used for a BENEFIT slider case it has to be taken that the lead wire stays absolutely straight along the intended direction of moving. Using the provided fixation screws the BENEFIT is provisionally tightened on the implantable screws and the accessories are connected and activated. Now the fixation screw is finally tightened (max. 10 Nm) on the implantable screws. When the both movement is occurring, the accessories have to be activated several times until the final tooth position is reached. The BENEFIT System used in tandem 2 screws connected with a BENEFIT slider is used in the palate. The anterior palate is the most suitable insertion area. The screws can be inserted medially or paramedial. The adaptation and connection of BENEFIT accessories should be done as soon as possible after the screw insertion either directly on the patient or in a laboratory using the impression cap (REF 33-54410) and the laboratory anvil (REF 33-54425) using a PVS impression.

**11. PILOT HOLE DRILLING**

Self-drilling miniscrews normally do not need a pilot hole prior to screw insertion in the bone with the following (but not limited to) exceptions. Insertion of 1.5 mm screws in the mandible requires a predrilling through the outer cortex. Palatally, a perforation of the outer cortical bone in adult patients is recommended. The diameter of the pilot hole depends on the screw diameter. The drills are colored and correspond to the color of the screw holders in the sterile blister. Green = 1.5 mm screw / 1.1 mm drill, red = 2.0 mm screw / 1.4 mm drill, grey = 2.3 mm / 1.8 mm drill. For each patient a new drill must be used. Ensure a maximum of 500 - 800 rpm and provide adequate cooling using a sterile, cooled, physiological saline solution. We recommend a manual drilling using a normal contra-angled hand piece with the manual turning unit 10-63025.

**12. ACCESSORIES**

The Quattro® and Quattro® Mini screws can be combined with commercially available orthodontic appliances like elastics, springs, wires, etc. The BENEFITs can be bend using commercially available bending pliers. General accessories for the BENEFIT System are shown in KAT-PB-032. Furthermore, refer to KAT-002 and KAT-003 for BENEFIT® and Quattro® System products respectively.

**13. DELIVERY CONDITIONS**

Implantable screws are delivered in sterile state. Sterilization is achieved by irradiation.

Re-sterilization of implantable screws is absolutely forbidden.

Do not use the sterile devices after expiration date.

Implantable screws must not be used if package has been opened or damaged. Device sterility may be jeopardized.

Abutments, fixation screw and plates, lingual sheath, mobilizers, springs and tubes are delivered in non-sterile state. They must be reprocessed before use.

Reuse of implantable screw and components, as well as the multiple reprocessing on components are absolutely forbidden.

**14. REPROCESSING PROCEDURES**

Recommendation: Components must be checked after reprocessing and before surgical procedure meticulously. Only cleaned components must be used.

Precaution: Pay attention to the manufacturer's detergent instructions regarding concentration, temperature and soaking time as well as post-rinsing. We assume the doctor has experience and knowledge of standard protocols regarding dental practice and limit these instructions only to the use of this product.

Precaution: Use freshly prepared solutions and demineralized water only.

PSM MEDICAL SOLUTIONS components are made out of suitable materials. There are no concerns regarding material resistance and/or any known sensitivity to process parameters during reprocessing (heat, cleaning agents etc.) which may affect the safety of our devices.

**15. AUTOMATED CLEANING AND DISINFECTION PROCEDURES**

STEP	DESCRIPTION
Preparation	Whenever possible disassemble the instruments. Place the disassembled instruments in the washer/disinfector (e.g. Miele & Cie. KG Gütersloh, Type: G 7336 CD). Ensure that the instruments do not touch each other. If products with narrow lumens or cavities cannot be connected, they must be placed in a way that water and cleaners can drain of completely. Then start the program (e.g., Shortened Vario TD Program)
Pre Cleaning	For 4 mins
Cleaning	For 10 mins at 55°C using an alkaline cleaner (e.g.: Neodisher MediClean forte [0,5%])
Neutralization	For 6 mins using a neutralizer (e.g. Neodisher Z [0,5%])
Final Rinsing	For 3 mins with deionized water
Thermal Disinfection	For 95 sec at 95°C
Remotion	After completion of the program, remove the instruments from the washer/disinfector
Inspection	Control and wrap the instruments immediately after removal from the washer/disinfector

**16. PACKAGING**

Sort the cleaned and disinfected instruments into the sterilization trays and package them in adequate sterilization packages (Single or double packaging) and / or sterilization containers meeting the following requirements:

- acc. to AAMI ISO 11607, suitable for steam sterilization (temperature resistance up to at least 137 °C (279 °F), sufficient vapor permeability), adequate protection of instruments or sterilization packaging from mechanical damage and regularly maintained according to the manufacturer's instructions (sterilization container).

**17. STERILIZATION**

Sterilize the components in accordance with parameters defined below:

Type	Process
Type	Steam sterilization performed in fractionated vacuum process (pre-vacuum process or gravity process with sufficient product drying)
Temperature	132 °C (270 °F) or 134 °C (273 °F) (plus tolerance acc. to ISO 17665)
Time	min. 5 mins (132°C) or min. 18 mins (134°C for prion inactivation)

**18. STORAGE**

PSM devices must be stored in a suitable environment. The storage room has to be dust-free, with low microbiological contamination, dark and free of temperature fluctuations.

Keep the devices away from sunlight

Keep the devices in dry area

**19. DISPOSAL**

The products are to be disposed of according to the local laws and regulations.

**20. QUALITY AND WARRANTY**

PSM ensures to the user an immaculate, inspected quality of all products. The responsibility for proper use remains with the user. National regulations and their restrictions have to be observed. All our products are designed and manufactured to the highest quality standards. PSM, as manufacturer of the products, excludes all warranty claims and takes no responsibility for direct or subsequent damage resulting from:

- use for the wrong purpose
- improper use, application or handling
- lack of training
- improper or repeated reprocessing or use
- combination with foreign products (BENEFIT)
- disregard of the instructions for use

**21. USE OF PSM ORIGINAL PRODUCTS**

BENEFITs and BENEFIT accessories are produced and designed to be used together with BENEFIT Screws and BENEFIT® accessories. No part of the system should be replaced with a different manufacturer's product, even if the product appears to be visually and dimensionally compatible or identical to the PSM product. The use of other manufacturers' products along with PSM products can involve incalculable risks and/or contamination of the material and misalignment of the implant to the instrument, thereby endangering the patient, user or third parties. For an optimized traceability and patient file documentation each package contains adhesive and detachable product data labels.

**22. USED SYMBOLS EXPLANATION**

	Manufacturer		Keep away from sunlight
	Consult the Instruction for Use		Keep dry
	English Description		Do not Re-Sterilize
	Symbol for "Airtic Number"		Do not Reuse
	Batch Code		Caution, consult accompanying documents
	Qty Quantity		Barcode for Reference Number
	Sterilized by Irradiation		Global Trade Item Number
	Do not use if packaging is open or damaged		US Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner
	Use by		European conformity to the essential requirements with notified body number
	Non-Sterile		

Fig. 2

**For further informations please contact:**

Manufacturer | Hersteller:  
psm MEDICAL SOLUTIONS  
Gewerbestraße 10 | 78594 Gunningen  
Telefon +49 (74 24) 9 75 15-0  
E-Mail info@psm.ms  
Web www.psm.ms

Distribution in USA:  
psm North America, Inc.  
8090 Weiskopf | La Quinta, CA 92253 USA  
Toll-free: 800-733-1622 | E-Mail service@psm-na.us  
Web www.psm-na.us

Notified body / Benannte Stelle:  
TÜV Nord Cert | www.tuev-nord.de  
CE0044